




## GDP-Mannose (GDP-Man) Kit



### Notes:

- All reagents and kit components should be stored at -20 °C until use
- This kit is intended for:
- GDP-Man kit containing substrates (GTP, Man) and enzyme (NahK15697, PFManC, PPA) is mini test kit *in situ* which is necessary step for large scale glycosylation.
- lower cost, continuous generation of up to 18 μmoles of nucleotide sugar *in situ* for use with sugar transferases (not included)
- conversion of sugar derivatives (not included) to the corresponding nucleotide sugar
- adding more GTP and Man to the reaction will continually generate large quantity of GDP-Man (not guaranteed) while NahK15697, PFManC and PPA are still active

## Quick start protocol

### Step 1: Inspect kit contents.

- Substrate Tube A: GTP (powder ; qty 1)
- Substrate Tube B: UTP (powder ; qty 1)
- Substrate Tube C: Mannose sugar (powder ; qty 1)
- Enzyme Tube D: NahK15697 (powder ; qty 1)
- Enzyme Tube E: PFManC (powder ; qty 1)
- Enzyme Tube F: EcPPA (powder ; qty 1)
- Reaction Tube G: Sterile empty tube (qty 1)
- Centrifuge all tubes briefly to pellet any material from walls of tube before opening tubes

For technical assistance: [sales@chemilyus.com](mailto:sales@chemilyus.com)



### **Step 2: Assemble additional components (not included).**

- Create a 20X Buffer Solution (1M Tris pH 8.0)      Add 300 $\mu$ L to Reaction Tube F.
- Create a 10X Salt Solution (200mM MgCl<sub>2</sub>)      Add 100  $\mu$ L to Reaction Tube F.

### **Step 3: Prepare reagents.**

- Add 100  $\mu$ L of dH<sub>2</sub>O from Step 2 to Substrate Tube A. Tap gently to mix. Centrifuge briefly to pellet any insoluble material. Transfer all 100  $\mu$ L to Reaction Tube G
- Repeat with Substrate Tube B
- Repeat with Substrate Tube C. [*Note: If using a Mannose derivative instead of Mannose as the sugar substrate, skip this step. Mannose derivative not included.*]
- Repeat with Enzyme Tube D [*Note: Enzymes should always be added to Reaction Tube G last*]
- Repeat with Enzyme Tube E
- Repeat with Enzyme Tube F

### **Step 4: Initiate nucleotide sugar reaction.**

- Ensure that Reaction Tube G contain all reagents. Final reaction volume is 1000  $\mu$ L
- Incubate Reaction Tube G for 3 h in 37 °C water bath
- After 3 h, centrifuge briefly to pellet condensation and any insoluble material
- Reaction Tube G now contains GDP-Mannose for glycosylation reactions

### **Step 5: Glycosyltransferase reaction.**

- Reaction Tube G can still actively produce GDP-Mannose *in situ*
- Add a Mannose transferase and target substrate to Reaction Tube G
- Incubate Reaction Tube G in 37 °C water bath for 24 hr to initiate glycosylation of the target substrate

[*Note: Glycosylation rate may vary by transferase and target substrate*]