




## UDP-glucuronic acid (UDP-GlcA) Kit



### Notes:

- All reagents and kit components should be stored at  $-20\text{ }^{\circ}\text{C}$  until use
- This kit is intended for:
- UDP-GlcA kit containing substrates (ATP, UTP, GlcA) and enzyme (AtGlcAK, AtUSP, PPA) is mini test kit *in situ* which is necessary step for large scale glycosylation.
- lower cost, continuous generation of up to 18  $\mu\text{moles}$  of nucleotide sugar *in situ* for use with sugar transferases (not included)
- conversion of sugar derivatives (not included) to the corresponding nucleotide sugar
- adding more ATP, UTP and GlcA to the reaction will continually generate large quantity of UDP-GlcA (not guaranteed) while AtGlcAK, AtUSP and PPA are still active

## Quick start protocol

### Step 1: Inspect kit contents.

- Substrate Tube A: ATP (powder ; qty 1)
- Substrate Tube B: UTP (powder ; qty 1)
- Substrate Tube C: GlcA sugar (powder ; qty 1)
- Enzyme Tube D: AtGlcAK (powder ; qty 1)
- Enzyme Tube E: AtUSP (powder ; qty 1)
- Enzyme Tube F: EcPPA (powder ; qty 1)
- Reaction Tube G: Sterile empty tube (qty 1)
- Centrifuge all tubes briefly to pellet any material from walls of tube before opening tubes.

For technical assistance: [sales@chemilyus.com](mailto:sales@chemilyus.com)



### **Step 2: Assemble additional components (not included).**

- Create a 20X Buffer Solution (1M Tris pH 8.0)      Add 300  $\mu$ L to Reaction Tube G.
- Create a 10X Salt Solution (200mM  $MgCl_2$ )      Add 100  $\mu$ L to Reaction Tube G.

### **Step 3: Prepare reagents.**

- Add 100  $\mu$ L of  $dH_2O$  from Step 2 to Substrate Tube A. Tap gently to mix. Centrifuge briefly to pellet any insoluble material. Transfer all 100  $\mu$ L to Reaction Tube G
- Repeat with Substrate Tube B
- Repeat with Substrate Tube C. *[Note: If using a GlcA derivative instead of GlcA as the sugar substrate, skip this step. GalNAc derivative not included.]*
- Repeat with Enzyme Tube D *[Note: Enzymes should always be added to Reaction Tube G last]*
- Repeat with Enzyme Tube E
- Repeat with Enzyme Tube F

### **Step 4: Initiate nucleotide sugar reaction.**

- Ensure that Reaction Tube G contain all reagents. Final reaction volume is 1000  $\mu$ L
- Incubate Reaction Tube G for 6 h in 37  $^{\circ}C$  water bath
- After 6 h, centrifuge briefly to pellet condensation and any insoluble material
- Reaction Tube G now contains UDP-GlcA for glycosylation reactions

### **Step 5: Glycosyltransferase Reaction.**

- Reaction Tube G can still actively produce UDP-GlcA
- Add a GlcA transferase (such as PmHAS ; Chemily product EN01025) and target substrate (such as GlcNAc) to Reaction Tube G
- Incubate Reaction Tube G in 37  $^{\circ}C$  water bath for 24 hr to initiate glycosylation of the target substrate

*[Note: Glycosylation rate may vary by transferase and target substrate]*



**Chemily Glycoscience**

For technical assistance: [sales@chemilyus.com](mailto:sales@chemilyus.com)